

# M800 vet

## Handheld Veterinary Monitor





SpO2  
+  
3-Lead  
ECG



### Parameter:

### Standard:

SpO2, Pulse Rate

### Optional:

3-Lead ECG

Nellcor SpO2

### Specifications:

\*Run time > 14h

\*2.4" color TFT display

\*Auto rotation function

\*Suitable for different types of animal; Apply to hospital clinic, vet shop

\*Easy to operate & carry convenient for clinical

// **Automatic power-off** function for power saving.

**Dual work modes:**  
continuous monitoring & spot check


# SPO2 -Clinical Application

- Helps establish the animals' **blood oxygen saturation**
- Helps view the pleth, which helps in understanding the **blood flow**
- Helps diagnose **hypoxia**
- Helps monitor the animals' **heart rate**



# ECG -Clinical Application

- An ECG allows an accurate determination of your animal's **heart rate**.
- It can also determine whether the heart is beating with a Normal pace, or whether it is beating too quickly or too slowly. An **elevated or decreased heart rate** can suggest specific medical issues that may warrant investigation.
- It can help diagnose **heart arrhythmias** that can be intervened and treated.
- It can help diagnose **ischemic changes in the heart**.
- It can help **diagnose electrolyte imbalance**.
- It can help diagnose **pericarditis and other heart conditions**
- 



# Recommendations for ECG monitoring



- ✓ Emergencies
- ✓ Anaesthesia
- ✓ Post surgical cases
- ✓ If an arrhythmia suspected on auscultation or history
- ✓ Drug toxicity
- ✓ Electrolyte disturbances
- ✓ To assess effectiveness of cardiac drugs

# M800 Vet Accessories



3-lead Snap Connector  
ECG Cable




ECG alligator clip Electrodes



BLT digital SpO2 ear clip  
reusable sensor +tongue clip



SpO2 extension cable



Nellcor Reusable  
dura-Y multi-site  
sensor



# M880Vet Accessories

(Optional)



LB-03 Lithium battery + charger stand



M800 Vet Carry bag



Silicone Case



# Thank You!

